


# DNA Barcoding

## *ELLS Troubleshooting Guide*




Teacher's  
Guide

Version 2.0

ELLS – European Learning Laboratory for the Life Sciences

# Experimental Troubleshooting

Version 2.0




# Experimental Troubleshooting, contd.

Version 2.0

<b>No band on gel</b>	See solutions for "No DNA"
<b>Faint band on gel</b>	See solutions for "Low DNA yield and/or quality"
<b>Poor sequencing quality</b>	<div> <div>Sample collection (all protocols)</div> <div> <div>Insufficient starting material</div> <div>Repeat DNA extraction using larger amount of starting material. Ensure your plant tissue is fresh (or has been stored at -20°C after picking it fresh). In the case of hard material (e.g. seeds), disrupt the plant tissue using pestle and mortar.</div> </div> </div> <div> <div>DNA extraction (Spin-column protocol only)</div> <div> <div>Insufficient binding of lysate to spin column membrane</div> <div>Ensure lysate volume prior to adding Buffer AW1 is accurately determined.</div> </div> <div> <div>Insufficient cell lysis</div> <div>Ensure sample is properly homogenised. Ensure sample is incubated sufficiently in the presence of lysis buffer AP1, at the appropriate temperature.</div> </div> <div> <div>DNA degraded during DNA extraction process</div> <div>Ensure the centrifugation step prior to loading lysate onto spin column is carried out.</div> </div> </div> <div> <div>DNA extraction (Direct PCR protocol only)</div> <div> <div>Insufficient cell lysis</div> <div>Ensure sample is properly homogenised. Ensure sample is sufficiently incubated in the presence of the extraction buffer.</div> </div> </div> <div> <div>PCR amplification (all protocols)</div> <div> <div>Insufficient amount of primer added to PCR reaction</div> <div>Repeat PCR using correct amount of forward and reverse amplification primer.</div> </div> <div> <div>Insufficient amount of DNA template added to PCR reaction</div> <div>Repeat PCR using correct amount of DNA template.</div> </div> <div> <div>Incorrect thermal cycling protocol used</div> <div>Repeat PCR using correct cycling protocol.</div> </div> <div> <div>Insufficient amount of PCR product loaded on gel</div> <div>Repeat gel electrophoresis loading larger amount of PCR product.</div> </div> </div> <div> <div>PCR amplification (Direct PCR protocol only)</div> <div> <div>No "Additives Mix" added to PCR reaction</div> <div>Repeat PCR using correct amount of "Additives Mix".</div> </div> </div> <div> <div>PCR clean-up (all protocols)</div> <div> <div>Insufficient clean-up, leading to PCR reagents interfering with sequencing reaction</div> <div>Repeat PCR clean-up</div> </div> </div>
<b>I have an issue which is not described here</b>	Please describe your issue online via the interactive Troubleshooting Guide or write us an email at <a href="mailto:emblog-ells@embl.org">emblog-ells@embl.org</a>


# Sequence Analysis Troubleshooting

Version 2.0

Sequence Analysis Troubleshooting 	
Issue with number of sequences returned	<p>I did not get any sequences back → No sequence analysis possible → For further solutions see "Experimental troubleshooting" &gt; "No DNA obtained" and "Poor sequencing quality"</p> <p>I only got one out of two sequences back → Use single seq. read for further analysis → For further solutions see "Experimental troubleshooting" &gt; "Poor sequencing quality"</p>
Issue with sequence coverage (length) and/or quality	<p>My sequence(s) are short and don't cover the whole barcode → Use short sequence(s) for further analysis → For further solutions see "Experimental troubleshooting" &gt; "Poor sequencing quality"</p> <p>My sequence(s) have poor quality (as seen in chromatogram) → Only use sequence(s) with reliable chromatogram(s) → For further solutions see "Experimental troubleshooting" &gt; "Poor sequencing quality"</p>
Issue with contig assembly	<p>Contig assembly is impossible → Use single seq. read for further analysis → For further solutions see "Experimental troubleshooting" &gt; "Poor sequencing quality"</p>
Issue with ENA search	<p>ENA search results are unclear or not as expected → See solutions at "ENA Search Troubleshooting"</p>
I have an issue which is not described here	<p>Please describe your issue online via the interactive Troubleshooting Guide or write us an email at <a href="mailto:emblog-ells@embl.org">emblog-ells@embl.org</a></p>

# ENA Search Troubleshooting

Version 2.0

ENA Search Troubleshooting 	
<b>Incorrect organism identification</b>	<p>Single species incorrectly identified</p> <ul style="list-style-type: none"> <li>Sequence information insufficient for database           <ul style="list-style-type: none"> <li>Ensure contig is assembled correctly</li> <li>To improve sequence coverage and quality see solutions to "Experimental Troubleshooting" &gt; "Poor sequence quality"</li> </ul> </li> <li>No ENA entry for organism           <ul style="list-style-type: none"> <li>No identification via ENA possible at the moment. Come back later to check for new entries. Use phenotypic characteristics of organism for species identification.</li> </ul> </li> </ul> <p>Expected species is only identified to genus level</p> <ul style="list-style-type: none"> <li>Sequence information insufficient for database           <ul style="list-style-type: none"> <li>Ensure contig is assembled correctly</li> <li>To improve sequence coverage and quality see solutions to "Experimental Troubleshooting" &gt; "Poor sequence quality"</li> </ul> </li> <li>No ENA entry for organism           <ul style="list-style-type: none"> <li>No identification via ENA possible at the moment. Come back later to check for new entries. Use phenotypic characteristics of organism for species identification.</li> </ul> </li> </ul>
<b>Multiple top hits (entries with identical ENA values)</b>	<p>Unable to distinguish between different species of the same genus</p> <ul style="list-style-type: none"> <li>Sequence information insufficient for database           <ul style="list-style-type: none"> <li>Ensure contig is assembled correctly</li> <li>To improve sequence coverage and quality see solutions to "Experimental Troubleshooting" &gt; "Poor sequence quality"</li> </ul> </li> <li>No ENA entry for organism           <ul style="list-style-type: none"> <li>No identification via ENA possible at the moment. Come back later to check for new entries. Use phenotypic characteristics of organism for species identification.</li> </ul> </li> </ul> <p>Unable to distinguish between different genera</p> <ul style="list-style-type: none"> <li>Sequence information insufficient for database           <ul style="list-style-type: none"> <li>Ensure contig is assembled correctly</li> <li>To improve sequence coverage and quality see solutions to "Experimental Troubleshooting" &gt; "Poor sequence quality"</li> </ul> </li> <li>No ENA entry for organism           <ul style="list-style-type: none"> <li>No identification via ENA possible at the moment. Come back later to check for new entries. Use phenotypic characteristics of organism for species identification.</li> </ul> </li> </ul>
<b>I have an issue which is not described here</b>	Please describe your issue via the interactive Troubleshooting Guide or write us an email at <a href="mailto:emblog-ells@embl.org">emblog-ells@embl.org</a>



