

Tryptophan fluorescence (WF) for total protein and peptide determination

Proteomics Core Facility
Tel +49 6221 387-8388
pcf@embl.de

EMBL Heidelberg
Meyerhofstrasse 1
69117 Heidelberg
Germany

www.embl.org

https://www.embl.de/proteomics/proteomics_services/index.html

Principle

The fluorescence spectrometry of tryptophan offers a simple, sensitive, and direct method for protein and peptide assays (Wiśniewski & Gaugaz, 2015). The WF assay is fully compatible with SDS and other solutes that are commonly used for the lysis of tissue and cells. The assay can be carried out on a standard fluorescence spectrometer with cuvettes and in a 96-well format using a plate reader. The method is particularly suitable for determination of peptide content in diluted samples.

Material

Reagents/Material	Details
8 M Urea in 10 mM HEPES pH 8.5	Prepare 48g of urea in 100 mL 10mM HEPES, pH 8.5 (NaOH) Note: Don't need to make it fresh
Albumin Standard (Pierce, A56979)	Pierce™ Dilution-Free™ BSA Protein Standards, multichannel pipette compatible. Ready to pipet, BSA: 0.125, 0.25, 0.5, 1, 2, 5, 10 mg/mL
96-well plate (Black, flat-bottomed, polystyrene plates)	For example: Corning CLS3915  Or similar plates Note: You can reuse the plates after washing them.

Sample Dilution

Sample	Dilution for assay
Total lysate or peptide	Dilute your sample in the assay buffer to end up in concentration range 0.125-10 mg/mL of the standard

Microplate Procedure

- Pipette 10 µL of each standard or your sample replicate into a **black** microplate well
 - Tip: Pipet your standard in the 1st, 2nd and 3rd column, then your samples rather than in rows because the standard has 8 concentrations (including 0 mg/mL)
- Add 200 µL of the urea assay buffer to each well and mix the plate thoroughly on a plate shaker for 30 seconds.
- Measure the fluorescence on a compatible instrument
Note: the signal is stable over hours

Settings

Important: Excitation at 295 nm and emission recorded between 320-400 nm and measurement at 20 °C

Here for example for a Tecan reader:

Mode	Fluorescence Reading	Top
Excitation Wavelength	295	nm
Emission Wavelength	355	nm
Excitation Bandwidth	5	nm
Emission Bandwidth	20	nm
Gain	100	Manual
Number of Flashes	100	
Flash Frequency	400	Hz
Integration Time	50	µs
Lag Time	0	µs
Settle Time	0	Ms
Z-Position (Manual)	20000	µm

Reference

- Wiśniewski JR, Gaugaz FZ. Fast and sensitive total protein and peptide assays for proteomic analysis. Analytical chemistry. 2015 Apr 21;87(8):4110-6. PMID: **25837572**. DOI: [10.1021/ac504689z](https://doi.org/10.1021/ac504689z)