

Laser	Channel Name-Collection Filter	Fluorochromes (Some)
355nM UV	355-379/28	BUV395 (Install the 390/22 BP instead)
	355-515/30	BUV496, -DAPI and Hoechst (Install the 450/50 BP instead)
405nM Violet	405-450/50	DAPI, Pacific Blue, AlexaFluor 405, BV421, V450, BFP, CFP, Cell Trace Violet
	405-525/50	Sytox Blue, Aqua Dead, AmCyan, V500, AlexaFluor 430, Pacific Orange, Qdot 525
	405-610/20	BV605, Qdot 605, eFluor 605NC, Fixable Yellow Dead, Pacific Orange
	405-660/20	BV650, Qdot 655
	405-710/50	BV711
	405-780/60	BV786
488nM Blue	488-530/30	FITC, GFP, YFP, AlexaFluor 488, mCitrine, Syto9, SYTOX Green, Calcein, Venus, DyeCycle Green
	488-710/50	PerCP-Cy5.5, PerCP, eFluor710, (DRAQ5)
561nM Yellow-Green	561-586/15	PE, DsRed, RFP, Tomato, SYTOX Orange
	561-610/20	PI, mCherry, Texas Red, PE-CF594
	561-670/14	7-AAD, PE-Cy5
	561-710/50	PE-Cy5.5., PE-Cy5, 7-AAD
	561-780/60	PE-Cy7, APC-Vio770, PE-Vio770
633nM Red	640-670/30	APC, AlexaFluor 647, DRAQ5, DRAQ7, Fixable Far Red, TO-PRO-3
	640-730/45	AlexaFluor 700, APC-Cy5.5, DRAQ7, DyeCycle Ruby
	640-780/60	APC-Cy7, AlexaFluor 750, APC-eFluor780, Fixable near IR, APC-H7